Department:	Dental	Number:	HHS-5105 D
		Attachment:	
Policy:	Operatory Disinfection		
Date: Februa	ry 1995 Revised: Ju	ly 2001, July 200	7, December 2012
Approved by:	Juig Ul		Date: 5/1/13

Policy Statement:

All treatment rooms must be cleaned and disinfected between patient use.

Purpose:

To eliminate any possible cross-contamination between patients and health care providers.

Procedure:

While wearing P.P.E's use Caviwipes towelette to completely preclean surface and remove all gross debris.

Spray Cavicide disinfectant spray directly onto surfaces

Allow to remain visibly wet for 30 seconds.

Use a second Caviwipe towelette to thoroughly wipe the wet surfaces leave wet for three minutes the time recommended by manufacturer.

Repeated use of the product may be required to ensure that the surface remains visibly wet for three minutes to be effective as a disinfectant.